

SILC  SKIN™

REMOVES WRINKLES WHILE YOU SLEEP

## SILCSKIN COLLETTE PADS

Instructions For Use



# { SILC SKIN }

Calvet Cosmetics, LLC | 14900 Magnolia Street | Suite 55186 | Sherman Oaks, CA 91413 | USA  
Phone: {855} 745.2756 | Fax: {818} 988.1346 | Email: info@silcskin.com | [silcskin.com](http://silcskin.com)

---



## COLLETTE PAD

*How the Collette Pad works and Instructions for Use.*

Our Collette Pads are composed of 100% Medical Grade Silicone. Electro static energy provided by the pads when worn is supposed to improve the irregular collagen structure and also creates an increased blood flow which results in an improved skin tone.

Silicone should contribute to moisturizing the skin by drawing the skin's moisture to the surface. This process has been shown to hydrate the skin and is an important factor in healing scars. We are simply treating wrinkles as scars by hydrating, improving the Collagen structure, and flattening them out. Prolonged use will help return neck and clavicle skin to a more youthful and supple appearance, by softening fine lines and wrinkles.

### INSTRUCTIONS FOR USE:

1. Make sure the neck and chest area is thoroughly clean, dry and free of moisturizers, sunscreens, etc. Sweating at night may interfere with adhesion. Dry well.
2. Remove pad from packaging. Remove the storage tray. Cut open poly bag containing the pad. Remove pad from backing.
3. The pad is self-adhesive. The pad will begin to lose adhesive qualities in approximately 30 days. Once this happens a new pad should be applied.
4. Apply the pad to the desired area of the neck or chest. You can put on the pad when going to bed, or any time of day. Use for an hour or more at a time.
5. Remove the pad slowly in the morning in an upwards motion. Place on the storage tray, and put the tray inside the packaging. Insert the tab on the package into the slit to close.
6. Pads should be washed occasionally to remove exfoliated skin cells that may interfere with the adhesive quality of the pad.
7. Pads should be washed with a soap that is free of moisturizers, emollients or oils. ie. Baby Shampoo. Allow pad to air dry face up on the storage tray in a clean location. Do not dry with a towel. Once dry, place face down on tray and return to package.
8. Some edge lifting is normal due to body contours. Do not use on sensitive, irritated, cut, scratched or burned skin. Do not use if there is a known sensitivity to silicone. Discontinue use if a rash or irritation occurs.
9. Our Collette pads are made of 100% medical grade silicone. There are no allergy causing additives or adhesives used. Pad is not transdermal.



Calvet Cosmetics, LLC  
14900 Magnolia Street  
Ste 55186, Sherman Oaks  
California 91413 USA

Phone: (855) 745.2756  
Fax: (818) 988.1346  
Email: [info@silcskin.com](mailto:info@silcskin.com)  
[silcskin.com](http://silcskin.com)

SILC  SKIN™  
REMOVES WRINKLES WHILE YOU SLEEP